

Coolshaping Comparison



Coolshaping vs Coolshaping2

Device		Coolshaping	Coolshaping2
Peltier dimension (W*H mm)	Large	217*64	186*76.5
	Medium	189*75	135.8*51.3
	Small	146*62	78.5*37
	mini	66*24	78.5*37
Handpiece size (D*W*H mm)	Large	217*125*141	225.5*120*115
	Medium	189*125*138	188.7*96*93
	Small	146*100*111	139*89*119
	mini	100*60*115	131*68*98
Cooling Temp.		-10 ~ 5°C	-10 ~ 5°C
Vacuum pressure		80~520 mmHg	80~520 mmHg
Vacuum pulse		4 pulse (CW, Low, Middle, High)	CW
Timer		1 ~ 60 min	1 ~ 60 min
Peltier location		Two sides	Whole internal part of the handpiece
Treatment time		45min	35min
Number of handpiece used simultaneously		2	4
Main unit size		380*443.5*1,067.2	498*774*1,119

Coolshaping vs Coolshaping2

(A)



(B)








(A) Coolshaping 1 device, (B) Coolshaping 2 Device




Coolshaping vs Coolshaping2

Components

(A)

Main Body	Medium Handpiece 2ea	Large COOL PAD 20ea (365 x 265)
		
Safety Key + Holder	Touch Pen	funnel hose
		
User's manual	Power cable	Fuse
		

Optional Components


Large Handpiece	Small Handpiece	Medium COOL PAD (260 x 230)	Mini Handpiece
			
Velcro for mini-handpiece attachment (4ea)	Mini COOL PAD (180 x 130)	Mini Handpiece liner	Mini handpiece Band (1Set)
			

(B)

Main Body	Large Handpiece 2	Medium Handpiece 1	Mini Handpiece 1
			
COOL PAD & LINER (L)	COOL PAD & LINER (M)	COOL PAD & LINER (MINI)	Bandage Set
			
COOL PAD : 20ea PAD LINER : 20ea	COOL PAD : 10ea PAD LINER : 10ea	COOL PAD : 10ea PAD LINER : 10ea	
Funnel Set	Body Cushion	Treatment Checking ruler Set	Touch Pen
			
Safety Key	Key Chain	The User's Manual	Power Cable
			
Option			
Small Handpiece			
			

(A) Coolshaping 1 device, (B) Coolshaping 2 Device


Coolshaping vs Coolshaping2

	L	M	S	Mini(Chin)
Coolshaping 1				
Coolshaping 2				


Coolshaping handpieces comparison

Coolshaping vs Coolshaping2

Coolshaping 1



Coolshaping 2



Coolshaping handpieces comparison